Gebruiksaanweisung

HIPP medical Plattensystem

Implantate: steril | unsteril

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HIPP medical Plating System

Implants: sterile | non-sterile

These Instructions for Use are not
For distribution in the USA.

Currently, not all products are
available in all markets.

© HIPP medical AG.
Subject to change.

HTG | **HIPP** medical AG
HIPP TECHNOLOGY GROUP

HIPP medical AG
Wilhelmstrasse 19
78600 Kolbingen
Germany

www.hipp-medical.com
info@hipp-medical.com
incident@hipp-medical.com

T: +49 7463 99 5167

CE 0483

Notified Body:
mdc medical device certification GmbH
www.mdc-ce.de
Identification number: 0483

Contents

Introduction Important Information for the User / Treating Physician / Third Parties.....	3
Material(s).....	3
Description.....	4
Intended Use / Clinical Benefits	4
Indications for the	5
Contraindications.....	5
Possible adverse effects (risks and side effects)	5
General adverse events	6
Notes and Warnings and Safety Instructions.....	7
Intraoperative phase	7
Postoperative Care	7
Compatibility and Instructions for Use	8
Disposable product / Single-use product	8
Packaging, delivery and storage.....	9
Disposal.....	9
Machine processing in the washer-disinfector (WD) and steam sterilization	10
Symbols for instructions for use and packaging/label	12
Procedure for serious incidents	12
Safety and Clinical Performance Summary (SSCP)	12
Warranty Statement.....	12

Marking of warnings and safety instructions in these operating instructions

Marking of notes in these operating instructions

Use and Care

Introduction | Important Information for the User / Treating Physician / Third Parties

! These instructions for use contain important warnings and precautions regarding the use of the **HIPP medical Plating System**. These instructions for use must be made available to the user / treating physician / third parties and the hospital staff at all times. The treating physician must be thoroughly familiar with the appropriate surgical technique. Sufficient experience in handling plates and screws is strongly recommended prior to using the **HIPP medical Plating System**.

! Read these instructions for use in their entirety and use the products only if you have understood these instructions.

⚠ Users, the treating physician, and third parties must be familiar with the surgical technique to be used. The treating physician is responsible for selecting the appropriate surgical technique for the specific patient and the procedure.

! **HIPP medical** cannot be held liable for complications arising from the use of the products that are beyond **HIPP medical's** control, including but not limited to product selection and deviations in application/handling and surgical technique.

⚠ The **HIPP medical Plating System** is intended to support the healing process and is NOT intended to bear weight or replace body structures in the event of an incomplete healing process.

Material(s)

The **HIPP medical Plating System** is manufactured from ELI-grade titanium alloy or unalloyed Grade 4 titanium:

Material(s)	Material Specification	Standard(s)
Titanium alloy	Ti-6Al-4V ELI	EN ISO 5832-3 / ASTM F136

Chemical composition according to ASTM F136:

Element	Mass fraction [%]
Aluminum (Al)	5.5 – 6.5
Carbon (C)	≤ 0.08
Iron (Fe)	≤ 0.25
Hydrogen (H)	≤ 0.012
Nitrogen (N)	≤ 0.05
Oxygen (O)	≤ 0.13
Titanium (Ti)	Residual
Vanadium (V)	3.5 – 4.5

Material(s)	Material Specification	Standard(s)
Unalloyed titanium	Titanium Grade 4	EN ISO 5832-2

Chemical composition according to EN ISO 5832-2:

Element	Mass fraction [%]
Carbon (C)	≤ 0.08
Iron (Fe)	≤ 0.50
Hydrogen (H)	≤ 0.0125
Nitrogen (N)	≤ 0.05
Oxygen (O)	≤ 0.40
Cobalt [Co]	≤ 0.10
Titanium (Ti)	Residual

Description

The **HIPP medical Plating System**, consisting of plates and screws, comprises implants for osteosynthesis and is used to fix bone fragments in cases of fractures, osteotomies, and arthrodeses. The **HIPP medical Plating System** is designed for use on the lower and upper extremities. It is available in sterile and non-sterile versions.

The **HIPP medical Plating System** ensures stable fixation of bone fragments, which promotes healing.

The use of these implants is not a substitute for natural tissue healing. The **HIPP medical Plating System** serves to further restrict the movement of a fractured bone and is intended only as an aid for fixing the fracture during the healing process.

The **HIPP medical Plating System**, consisting of plates and screws, includes the following models:

Fracture and Fusion Plates	Fracture Plates  Fusion plates 
Osteotomy plates	
Locking Cortical Screw	
	MVA screw 
	Screw with lentil head 
Cancellous screw	Partially threaded screw 
	Fully threaded screw 

Intended Use / Clinical Benefits

The **HIPP medical Plating System** is intended for the fixation of fractures, osteotomies, and fusions of the upper and lower extremities.

The **HIPP medical Plating System** consists of surgically invasive implants for the stabilization of bone fragments. They are intended exclusively for use by qualified medical professionals. The products are single-use devices for long-term use (>30 days) and should be removed once healing is complete.

Indications for the

The **HIPP medical Plating System** is indicated for:

- Fixation of bone fractures and bone fragments
- Stabilization of bones in arthrodesis procedures
- Fixation of bones in corrective and reconstructive osteotomies
- Bone reconstruction procedures

Contraindications

The selection of the appropriate implant must take the patient's general condition into account. The following factors may prevent or reduce the likelihood of a successful procedure; the **HIPP medical Plating System** should not be used in the following situations:

- Situations in which internal fixation is contraindicated for other reasons or in patients with health impairments (e.g., inflammation, immune defects, system diseases and metabolic disorders, insufficient bone quantity or quality, missing bone substance, circulatory disorders, increased leukocyte count) and/or lack of cooperation (e.g., alcoholism).
- Severe muscular, neurological, or vascular deficits in the affected limb.
- Active, acute, latent, potential, or chronic infections.
- All concomitant diseases that may endanger/impair the fixation and/or the success of the procedure or may increase the occurrence of postoperative complications.
- Fracture near the surgical area.
- Drug addiction, nicotine consumption, which can endanger the success of the procedure/surgery due to delayed bone/wound healing.
- Obesity, significant obesity of a patient can lead to additional stress and deformation of the product, which could accelerate fatigue of the material from which the implant is made and lead to deformation or destruction of it.
- Patients with conditions (e.g., mental illness) that limit their ability or willingness to follow postoperative instructions.
- Patient refusing to follow postoperative recommendations; advanced age or addiction (these conditions may cause the patient to ignore the restrictions and precautions during use of the implants).
- Large physical activities associated with strong vibrations that can lead to overloading of the implants and/or patient with high activity levels.
- If there is a proven intolerance/allergy to the material used
- Availability of conservative treatment options.
- Surgical procedures other than those listed in the indications.

Possible adverse effects (risks and side effects)

- Bone fragmentation or the so-called stress-shielding effect, which causes bone loss above, below or at the surgical site.
- Inadequate bone healing (delayed healing/non-healing/malunion)
 - Delayed union: bone fracture that takes longer than usual to heal.
 - Non-union: fractured bone fails to heal properly.
 - Malunion: fractured bone heals in an incorrect position (e.g., pseudoarthrosis, deformity)
- Bone fracture / refracture
 - Iatrogenic fracture: bone fracture that occurs as a direct result of the medical procedure itself.
 - Refracture: occurrence of a new fracture at the site of a previously healed fracture or a fracture in healing process
- Neurovascular impairment (nerve injury/nerve lesion/vascular injury/ischemia)
 - Nerve injury: Damage to nearby nerves during the insertion of the implants, leading to symptoms like numbness, tingling, or loss of motor function.
 - Nerve lesion: A nerve lesion in implants treatment refers to damage or injury to a nerve that occurs during the insertion, positioning, or removal of implants.
 - Vascular injury: Damage to blood vessels, which can result in bleeding, hematoma formation, or compromised blood flow to the affected area.

- Ischemia: Reduced blood supply to tissues, which can cause pain, pallor, and in severe cases, tissue death.
- Soft tissue injury (tendon rupture)
 - Soft tissue injury refers to damage to the muscles, tendons, ligaments, or skin that can occur during the insertion, positioning, or removal of osteosynthesis implants.
- Pain, discomfort, abnormal sensations, or palpability due to the presence of the product
 - Manipulation of a fracture is painful, so this is usually carried out using local, regional, or general anesthesia
 - Post-operative pain: Can occur e.g., from infections, soft tissue irritation. Effective pain management strategies are crucial to help patients recover comfortably.
 - Increased fibrous tissue reaction around the fracture site and/or the implant
- Postoperative wound complications
 - Improper wound healing: Factors that can impair wound healing and lead to complications include bacterial infection, necrotic tissue, foreign bodies, diabetes, smoking, malignancy, malnutrition, poor blood supply, global hypotension, hypothermia, immunosuppression (including steroids), emergency surgery, ascites, severe cardiopulmonary disease, and intraoperative contamination.
- Infectious diseases (osteomyelitis)
 - An infection is the invasion and multiplication of harmful microorganisms, such as bacteria, viruses, or parasites, within a host organism. Infection may be superficial or deep, which could lead to failure of fixation, malunion and osteomyelitis.
- Avascular necrosis (osteonecrosis)
 - Avascular necrosis (AVN), also known as osteonecrosis, is a condition where bone tissue dies due to a lack of blood supply.
- Displacement/instability of hardware (migration/exposure of metal parts/implant loosening)
 - Migration: movement of the implant from its original position.
 - Exposure of metal ware: implant becomes visible through the skin.
 - Loosening: implant becomes unstable and loses its secure position within the bone
- Implant failure (fractures/notches and pressure marks)
 - Breakage during and after treatment.
 - Repeated bending back and forth will cause fatigue or breakage of the implant.
 - Notches and pressure points also significantly reduce the mechanical strength of the implant
- Joint immobility
 - Inability to move a joint normally due to stiffness or restriction.

In general, the good material compatibility of implants made of titanium alloy Ti-6Al-4V ELI is supported by clinical data. Nevertheless, the following complications may occur:

- Fragment dislocations due to incorrect indication
- Neurovascular injuries due to surgical trauma
- Inflammatory reactions/hypersensitivity reactions
 - Inflammation is the body's natural response to injury or infection. Inflammation can occur without infection, such as in response to surgical trauma, but infection almost always causes inflammation.
 - A hypersensitivity reaction is e.g., hardware intolerance.
 - Allergic reactions
 - Foreign body reactions
- Infections can lead to failure of the procedure

General adverse events

In addition to the side effects, there is always the possibility of complications associated with any surgical procedure that are not necessarily related to the implant.

In addition to the failure of surgical procedures, many possible reactions can occur, but some of the most common are: Complications resulting from anesthesia and patient positioning (e.g., nausea, vomiting, neurological impairments, etc.), thrombosis, embolism, infection or injury of other critical structures, including blood vessels, excessive bleeding, damage to soft tissue, including swelling, cosmetic complications such as abnormal scar formation/skin necrosis, functional impairment of the musculoskeletal system, nerve injuries, postoperative wound complications and postoperative pain, discomfort or abnormal sensations due to the device, allergy or hypersensitivity reactions, loosening, bending or fracture of the device, healing in malposition, pseudarthrosis or delayed fracture healing, which may lead to implant failure and a repeat surgery.

Notes and Warnings and Safety Instructions

In addition to the notes highlighted in the individual chapters of these instructions for use, the following applies:

- ! The **HIPP medical Plating System** should only be accepted if the packaging and labeling are undamaged.
- ! Do not use the **HIPP medical Plating System** in procedures requiring a permanent implant.
- ! The attending physician should carefully assess bone quality before performing orthopedic procedures on patients.
- ! The **HIPP medical Plating System** should not be used if damage has been detected prior to and/or during use or implantation.
- ! The factors described in "Indications" and "Contraindications" should be taken into account when selecting the patient.
- ! The selection of the shape and size of the implant for each individual patient is a critical element that ensures the success of the surgery. The physician responsible for the selection is the attending physician.
- ! Preoperative and postoperative procedures, including knowledge of surgical techniques, along with correct positioning of the implants, are important factors that determine the success of the surgery.

Intraoperative phase

- ! Improper selection, placement, positioning and fixation of the implant may result in an undesirable result later.
- ! The **HIPP medical Plating System** must not come into contact with objects that could damage its surface. The plates and screws must not be mechanically processed or otherwise modified.
- ! Plates and screws of the **HIPP medical Plating System** that have come into contact with bone tissue, soft tissue, blood, synovial fluid and skin must not be reused and must be disposed of by clinic staff in accordance with the guidelines for the disposal of contaminated products. Residues of contamination on the implants can cause injury or infection to the patient, user or third parties.
- ! The mechanical processing of the **HIPP medical Plating System** (plates and screws) by the user is not permitted.
- ! It is essential to avoid deformation of the **HIPP medical Plating System** (plates and screws). Such tampering, in conjunction with improper use, can lead to surface defects and/or structural changes in the material, ultimately resulting in product failure.
- ! The use of **HIPP medical Plating System** in an MRI environment carries various risks unless the implant is marked as "MR safe" or "MR-related." These include, but are not limited to, heat generation and/or migration of the implant, and artifacts caused by the implant. However, under all circumstances, the user / attending physician should familiarize himself or herself with the contraindications and warnings of the manufacturer of the MRI device with which the examination is planned to be performed.
HIPP medical Plating System are not tested for MRI compatibility.

Postoperative Care

- ! The attending physician should give the patient instructions on postoperative behavior and gradual weight-bearing progression. Excessive load and/or activity can cause the **HIPP medical Plating System** (plates and screws) to break or bend, which may result in failure of the plates and/or screws and/or the treatment.
- ! The final decision regarding the timing of removal of the **HIPP medical Plating System** is the responsibility of the attending physician. It is recommended - if possible and applicable for the individual

patient - to remove the fixation products after the healing process is fully completed. This is especially true with young and active patients. The risk of adverse effects such as secondary infection, allergies, material fatigue fractures, implant failure and/or impaired blood circulation increases with the duration of the implant in the body.

! The patient should be tested for infections at regular intervals for the duration of the implant in the body.

! The **HIPP medical Plating System** is designed for temporary use and should be removed after healing.

Compatibility and Instructions for Use

The plates and screws of the **HIPP medical Plating System** are intended for use together and have been designed, validated, and approved exclusively as a complete, dedicated system. Combination with components outside this system is not permitted; only the plates and screws belonging to the **HIPP medical Plating System** may be used together.

! The use of the special applicator is mandatory.

The following tools are required for proper use:

- appropriate screwdriver / bit

Screw	Bit	Color	Diameter
Cancellous	Torx 10	Gold	Ø 4 mm
		Gray	Ø 4 mm
Cortical	Torx 10	Gold	Ø 3.5 mm
		Red	Ø 3 mm
	Torx 8	Blue	Ø 2.8 mm
		Green	Ø 2.5 mm
Locking cortical	Torx 10	Gold/Orange	Ø 3.5 mm
		Purple	Ø 4 mm
		Red	Ø 3 mm
	Torx 6	Gold	Ø 2 mm
	Torx 8	Blue	Ø 2.8 mm
		Green	Ø 2.5 mm

Disposable product / Single-use product

The plates and screws of the **HIPP medical Plating System** are single-use products.

Do not re-use (labeled on the package)

! Products intended for single use must not be reused. Reuse or processing (e.g., cleaning and re-sterilization) of disposable products may impair the structural integrity of the product and/or cause product failure, resulting in injury, illness, or death of the patient. Furthermore, the reuse or clinical processing of single-use products increases the risk of contamination, for example through the transmission of pathogens from patient to patient. This may also result in injury, illness or death of the patient, user or third party.

! Sterile **HIPP medical Plating System** can be reprocessed according to the following specification if they become non-sterile during surgery and have not yet been modified, have not yet been implanted and have not come into contact with bone tissue, soft tissue, blood, synovial fluid and skin.

! **HIPP medical** does not recommend the clinical processing of contaminated **HIPP medical Plating System**. **HIPP medical** products contaminated with bone tissue, soft tissue, blood, synovial fluid, and skin should not be reused under any circumstances and should be disposed of in accordance with the clinic's guidelines and rules. Even if components appear externally intact after use, minor defects and invisible material damage can cause material fatigue.

! Plates and screws from the **HIPP medical Plating System** that have already implanted must not be reused since the risks of material fatigue damage and infection are very high. **HIPP medical** accepts no liability in the event of non-compliance.

Packaging, delivery and storage

The plates and screws of the **HIPP medical Plating System** are supplied both sterile and non-sterile. The implants (plates and screws) are individually packaged and delivered in protective packaging. Upon delivery, the packaging should be undamaged.

Product is supplied sterile, sterilized by irradiation (marking on label)

Non-sterile / non-sterile product (marking on label)

- ! Store sterile plates and screws of the **HIPP medical Plating System** in their sterile packaging in a closed area, protected from sunlight, dust, moisture, and temperature fluctuations. Storage temperature range: +15°C to +30°C.

There is no limit to the storage life of non-sterile products. The products are made of non-degradable material, so there is no risk of loss of stability when stored under the recommended conditions. Store non-sterile plates and screws of the **HIPP medical Plating System** in their primary packaging in a closed area, protected from sunlight, dust, and moisture.

- ! Always check the expiration date before using sterile products. Do not use sterile products past their expiration date.

Can be used until / expiration date (marking on label)

- ! Check the integrity of the sterile packaging before use. Plates and screws from the **HIPP medical Plating System** should only be accepted if the packaging and labeling are undamaged.

Do not use if the packaging is damaged

Non-sterile **HIPP medical Plating System** must be cleaned, disinfected, and sterilized prior to surgical use. Prior to cleaning, the original packaging is to be completely removed. Prior to steam sterilization, pack the products in a validated packaging system (sterilization wrap or sterilization container).

In addition to the label, an auxiliary label may be affixed to the outer packaging, which contains specific requirements of a particular market sector (e.g., legal requirements of the country in which the product is to be marketed).

Disposal

The plates and screws of the **HIPP medical Plating System** are made of unalloyed Grade 4 titanium or Grade 4 titanium alloy ELI-grade.

HIPP medical Plating System contaminated with blood, tissue and/or body fluids and substances or compression screws that are no longer usable must not be reused under any circumstances and must be handled in accordance with the clinic's guidelines and rules of operator. Please ensure that the contaminated **HIPP medical Plating System** has been properly decontaminated (cleaned and sterilized) prior to disposal. If this is not done, they must be disposed of as infectious special waste.

Contaminated or unusable medical devices must be disposed of as (clinical) healthcare waste in a safe manner and in accordance with local healthcare guidelines and government and regulatory laws or guidelines.

When separating, recycling, or disposing of packaging materials, local government and regulatory laws on packaging and packaging waste must be observed where applicable.

The plates and screws of the **HIPP medical Plating System** are individually packaged with:

Specifications	Material
NON-STERILE	
Tip protector	TPE Thermoplastic Elastomer
Primary packaging	Polyamides/Polyethylene
Outer packaging	PE Polyethylene
Instructions for use	Thin printing paper
STERILE	
Tip protector	TPE Thermoplastic elastomer
Primary packaging	Polyamides/Polyethylene
Secondary packaging	Polyamides/Polyethylene
Outer packaging material / cardboard – non-sterile	100% primary wood fibers
Instructions for use	Thin printing paper

Machine processing in the washer-disinfector (WD) and steam sterilization

Validated procedure	Thermal disinfectant and steam sterilization
WARNINGS	Manual cleaning is not intended. For thermostable medical devices, machine processing with thermal disinfection followed by steam sterilization should be preferred over other procedures. Processing should only be carried out by trained personnel.
Recommendation and limitation of reprocessing	The cleaning disinfectors for processing must comply with national regulations (e.g., DGHM or FDA approval or CE marking, EN ISO 15883). For avoiding cross-contamination, assembly and disassembly should only be carried out by trained personnel with clean gloves. In the event of visual and/or functional damage, the plates and screws of the HIPP medical Plating System must not be reused. The product lifespan of plates and screws of the HIPP medical Plating System is limited to three (3) reprocessing cycles.
Identification and traceability	Items REF HM-Mxxx do not have direct marking; all other items in the product group are directly marked as follows: Product marking: Laser marked UDI Data Matrix with GTIN Batch number: Laser marked UDI Data Matrix with LOT CE marking: CE 0483 laser marked Fully threaded screw are not directed marked with CE due to lack of space. Traceability: Laser marked UDI Data Matrix
Accessories	Plates and screws of the HIPP medical Plating System must be placed in a suitable sieve tray. This tray is not included in the standard delivery.

INSTRUCTIONS FOR PROCESSING

Delivery	Non-sterile plates and screws of the HIPP medical Plating System must be processed prior to first use. The packaging of the implants must be carefully and completely removed prior to cleaning, sterilizing, and surgery. Never sterilize in the sterilizer with the packaging.
Place of use	For avoiding cross-contamination, assembly and disassembly should only be carried out by trained personnel with clean gloves .
Storage and transport	Plates and screws of the HIPP medical Plating System must be placed in a suitable sieve tray. A prerequisite for effective machine processing is the correct loading of sieve trays, inserts, holders, etc. in accordance with rinsing requirements. The sieve trays must not be overloaded, so that the products are well rinsed. Large-area products must be placed on the sieve tray in such a manner that they do not obstruct the cleaning of other products due to "shadowing." The products must be deposited or stored according to their mechanical sensitivity in such a manner that damage can be ruled out.
Pre-cleaning	If plates and screws of the HIPP medical Plating System are visibly contaminated, remove them with a lint-free damp cloth . No further pre-cleaning is required.

Preparation for decontamination	Make sure that t plates and screws of the HIPP medical Plating System are securely seated in the sieve basket.																												
Cleaning, mechanical / disinfection	During the validation of the automated cleaning process using washer-disinfectors (WD), the following minimum cycle parameters were applied: <table><tr><th>Step</th><th>Temperature [°C]</th><th>Time [min]</th><th>Cleaning agent</th></tr><tr><td>Pre-flushing</td><td>City water (cold)</td><td>1</td><td>-</td></tr><tr><td>Cleaning</td><td>55 ± 5</td><td>5</td><td>Alkaline cleaner Neodisher® MediClean forte 0.5%</td></tr><tr><td>Neutralization</td><td>DI water (cold)</td><td>2</td><td>Acid-based neutralizing agent Neodisher® Z 0.1%</td></tr><tr><td>Intermediate rinse</td><td>DI water (cold)</td><td>1</td><td>-</td></tr><tr><td>Thermal disinfection</td><td>90 (A0 3000)</td><td>5</td><td>-</td></tr><tr><td>Drying</td><td>60 ± 5</td><td>30</td><td>-</td></tr></table> Only WDs that comply with the general requirements for washer-disinfectors (WDs, described in part 1 of EN ISO 15883) should be used. In order to achieve an optimum cleaning result, process parameters such as dosing quantity and dosing temperature must be determined according to the specifications of the manufacturer of the cleaning agent. There are a number of points that require special attention during mechanical processing: • A prerequisite for effective machine processing is the correct loading of sieve trays, inserts, holders, etc. in accordance with rinsing requirements. • The sieve trays must not be overloaded, so that the instruments are well rinsed. • Large-area instruments must be placed on the sieve tray in such a manner that they do not obstruct the cleaning of other instruments due to "shadowing". • The instruments must be deposited or stored according to their mechanical sensitivity in such a manner that damage can be ruled out. • Products with cavities must also be completely rinsed on the inside. Special inserts with rinsing devices adapted to the products must be used for this purpose. • The items to be washed must be removed from the machine immediately after the program has ended, as remaining in the closed machine can cause corrosion due to residual moisture. • If elevated chloride concentrations are present in the water, pitting and stress corrosion cracking may occur on the instruments. Such corrosion can be minimized by using alkaline cleaners or demineralized water.	Step	Temperature [°C]	Time [min]	Cleaning agent	Pre-flushing	City water (cold)	1	-	Cleaning	55 ± 5	5	Alkaline cleaner Neodisher® MediClean forte 0.5%	Neutralization	DI water (cold)	2	Acid-based neutralizing agent Neodisher® Z 0.1%	Intermediate rinse	DI water (cold)	1	-	Thermal disinfection	90 (A0 3000)	5	-	Drying	60 ± 5	30	-
Step	Temperature [°C]	Time [min]	Cleaning agent																										
Pre-flushing	City water (cold)	1	-																										
Cleaning	55 ± 5	5	Alkaline cleaner Neodisher® MediClean forte 0.5%																										
Neutralization	DI water (cold)	2	Acid-based neutralizing agent Neodisher® Z 0.1%																										
Intermediate rinse	DI water (cold)	1	-																										
Thermal disinfection	90 (A0 3000)	5	-																										
Drying	60 ± 5	30	-																										
Visual inspection	After cleaning, the products must be visually checked for cleanliness. Products that have not become clean must be subjected to renewed cleaning (if necessary, manual pre-cleaning).																												
Preparation and packaging	After cleaning and disinfection, the products must be placed in packaging and trays suitable for sterilization in accordance with EN ISO 11607.																												
(Steam) sterilization	Only steam is permitted as a sterilization method. This must be able to come into contact with all points of the products in order to ensure complete sterilization. Sterilize the products according to the general hospital procedures for steam sterilization. Pack each product in packaging suitable for sterilization. Sterilization containers were used for sterilization validation in accordance with DIN EN ISO 11607-1 and DIN EN 868-8. The information (instructions for use) provided by the sterilizer manufacturer must be observed. Steam sterilization, fractional vacuum <table><tr><td>Temperature:</td><td>Min. 134°C.</td></tr><tr><td>Holding time:</td><td>Min. 5 minutes</td></tr><tr><td>Fractionation / pre-vacuum</td><td>Yes (minimum 3 pre-vacuum stages)</td></tr><tr><td>Drying:</td><td>Min. 10 minutes</td></tr></table> Allow products to cool to room temperature. Plates and screws of the HIPP medical Plating System that have not come into contact with blood, soft tissue, bones or body fluids may be sterilized and used again after machine processing. Plates and screws of the HIPP medical Plating System may be sterilized up to three (3) times.	Temperature:	Min. 134°C.	Holding time:	Min. 5 minutes	Fractionation / pre-vacuum	Yes (minimum 3 pre-vacuum stages)	Drying:	Min. 10 minutes																				
Temperature:	Min. 134°C.																												
Holding time:	Min. 5 minutes																												
Fractionation / pre-vacuum	Yes (minimum 3 pre-vacuum stages)																												
Drying:	Min. 10 minutes																												
Storage	Plates and screws of the HIPP medical Plating System are to be stored in a dry, clean, and dust-free environment. Protect from direct sunlight. Store sterile instruments in the sterile packaging in a closed cabinet, protected from dust, moisture, and temperature fluctuations. Temperature range for storage: +15°C to +30°C. The storage period of non-sterile products is not limited. The products are made of non-degradable material, such that there is no loss of stability due to storage under the recommended conditions.																												

Additional information

When processing multiple instruments in a reprocessing cycle, the maximum load of the particular WD or sterilizer must not be exceeded.

The preceding instructions have been validated as 'SUITABLE' by the medical device manufacturer for preparing a medical device for reuse. It is the responsibility of the reprocessor to ensure that the actual processing with the equipment, materials and personnel used in the processing facility achieves the desired results. State-of-the-art technology and often national laws require these processes and procedures to be validated and regularly monitored.

Neodisher is a registered trademark of Chemische Fabrik Dr. Weigert GmbH & Co. KG

Symbols for instructions for use and packaging/label

Symbol	Explanation	Symbol	Explanation
	Catalogue number		Unique product identification (UDI code)
	Isolated quantity		Manufacturer
	Non-sterile (Product is delivered non-sterile)		Do not use if the packaging is damaged and follow the instructions for use
	Do not reuse (single use only)		Store in a dry state
	Do not re-sterilize		Protect from sunlight
	Attention		Product is delivered sterile, radiation sterilized
	Observe instructions for use		Usable until
	Medical device		Temperature limits (+15°C to +30°C)
	CE marking and identification number of the notified body		Simple sterile barrier system with internal protective packaging
	Batch designation		Authorized representative in Switzerland ¹
	Date of manufacture		

¹ Applicable only to certain products in the product group

Procedure for serious incidents

All serious incidents related to **HIPP medical Plating Systems** must be reported to **HIPP medical** (incident@hipp-medical.com) and the competent national authority.

Safety and Clinical Performance Summary (SSCP)

The product group has the Basic UDI-DI **42556531Plates+Screws6W**. A Summary of Safety and Clinical Performance (SSCP) in accordance with Regulation (EU) 2017/745, Article 32, can be viewed using the Basic UDI-DI at the following web address: <http://ec.europa.eu/tools/eudamed>.

Warranty Statement

HIPP medical AG is solely responsible for ensuring that each individual product has been manufactured, inspected, and packaged with the utmost care. Since **HIPP medical** AG has no influence over or control of indications and/or applications, **HIPP medical** AG cannot be held responsible for complications or the failure of an application. **HIPP medical** AG individual products and sets are compatible with one another. Nevertheless, the user is required to ensure the compatibility of the products with one another prior to use. This applies particularly when the user employs **HIPP medical** AG products in conjunction with products from other manufacturers. Employees of **HIPP medical** AG are not authorized to modify the aforementioned conditions, extend liability, or assume additional product-related obligations.